

QuEChERS



Quick, Easy, Cheap, Effective, Rugged and Safe!
dSPE Tubes and standards for pesticide residue analysis from food products

See page 10 for details.



Buy one 10-pk. of Dual Vespel® Ring Inlet Seals at 25% off & Receive a 2nd 10-pk. FREE!



See page 8 for details.

Subscribe online!
www.restek.com/essentials

What's New!

PLOT Columns	4
Merlin Microseal Septa Kits	9
Canister Cleaning System	13
Ultra II LC Columns	14

Great Savings!

Rxi® GC Columns	2
Certified Autosampler Vial Kits	5
Generators	7
FREE Base Plate	7
Inlet Seals	8
Extraction Cells for ASE®	11
HPLC Detector Lamps	15
Leak Detector	16



Rxi[®] 3-in-1 Technology

Highest Inertness • Lowest Bleed • Exceptional Reproducibility


Exceptional Reproducibility

Engineered to assure column-to-column reproducibility.



Rxi[®]-5ms

- | | |
|---------------------|-------------------|
| 1. 1,6-hexanediol | 5. tridecane |
| 2. 4-chlorophenol | 6. 1-undecanol |
| 3. methyl nonanoate | 7. acenaphthylene |
| 4. 1-decylamine | 8. pentadecane |



Column: Rxi[®]-5ms, 30m, 0.25mm ID, 0.25 μ m (cat.# 13423) For conditions, visit www.restek.com and search for GC_EV00819

GC_EV00819

Make the Switch to Rxi[®] columns!

Rxi[®] Equivalency

	Rxi [®] -5ms	Rxi [®] -5Sil MS	Rxi [®] -1ms
Restek	Rtx-5	Rtx-5Sil MS	
Agilent	HP-5, HP-5ms, DB-5, Ultra-2	DB-5ms	HP-1, HP-1ms, DB-1, DB-1ms
Varian	CP-Sil 8 CB	VF-5ms, CP-Sil 8 CB LB/MS	VF-1ms, CP-Sil 5 CB LB/MS
Alltech	AT-5	AT-5ms	AT-1
Phenomenex	ZB-5	ZB-5ms	ZB-1, ZB-1MS



phases available

- Rxi[®]-1ms
- Rxi[®]-5ms
- Rxi[®]-5Sil MS
- Rxi[®] guard/retention gap columns
- Rxi[®]-35Sil MS
- Rxi[®]-17
- Rxi[®]-XLB

Rxi[®]-1ms Columns (fused silica)

(Crossbond[®] 100% dimethyl polysiloxane) Temp. range: -60°C to 330/350°C

Description	cat.#
0.25mm ID, 0.25 μ m, 30-Meter	13323-R25
0.25mm ID, 0.25 μ m, 60-Meter	13326-R25
0.32mm ID, 0.25 μ m, 30-Meter	13324-R25
0.32mm ID, 1.00 μ m, 30-Meter	13354-R25
0.32mm ID, 1.00 μ m, 60-Meter	13357-R25

Rxi[®]-5ms Columns (fused silica)

(Crossbond[®] 5% diphenyl/95% dimethyl polysiloxane) Temp. range: -60°C to 330/350°C

Description	cat.#
0.25mm ID, 0.25 μ m, 30-Meter	13423-R25
0.25mm ID, 0.25 μ m, 60-Meter	13426-R25
0.25mm ID, 0.50 μ m, 30-Meter	13438-R25
0.25mm ID, 1.00 μ m, 30-Meter	13453-R25
0.32mm ID, 0.25 μ m, 30-Meter	13424-R25

Rxi[®]-5Sil MS Columns (fused silica)

(Crossbond[®], selectivity close to 5% diphenyl/95% dimethyl polysiloxane) Temp. range: -60°C to 350°C

Description	cat.#
0.25mm ID, 0.25 μ m, 30-Meter	13623-R25
0.25mm ID, 0.25 μ m, 60-Meter	13626-R25
0.25mm ID, 0.50 μ m, 30-Meter	13638-R25
0.25mm ID, 1.00 μ m, 60-Meter	13697-R25
0.32mm ID, 0.50 μ m, 30-Meter	13639-R25

Lowest Bleed & Most Stable Bonded Carbowax® Column Available!


Stabilwax® Columns

(Crossbond® Carbowax® polyethylene glycol)

- Excellent inertness for active compounds.
- 260°C maximum operating temperature.

ID	df (μm)	temp. limits	length	cat. #
0.25mm	0.25	40 to 250/260°C	30-Meter	10623
0.32mm	0.25	40 to 250/260°C	30-Meter	10624
0.32mm	0.50	40 to 250/260°C	30-Meter	10639
0.32mm	1.00	40 to 240/260°C	30-Meter	10654
0.53mm	1.00	40 to 240/260°C	30-Meter	10655

Stabilwax® Structure



manufacturing
procedure
Guaranteed column-to-column
reproducibility

Unique Selectivity for Halogenated Compounds


Rtx®-200 Columns

(Crossbond® trifluoropropylmethyl polysiloxane)

- Temperature range: -20°C to 340°C
- Equivalent to USP G6 phase
- Similar phases: DB-200, DB-210

ID	df (μm)	temp. limits	length	cat. #
0.25mm	0.25	-20 to 320/340°C	30-Meter	15023
0.25mm	1.00	-20 to 290/310°C	30-Meter	15053
0.32mm	0.50	-20 to 310/330°C	30-Meter	15039
0.32mm	1.00	-20 to 290/310°C	30-Meter	15054
0.53mm	1.00	-20 to 290/310°C	30-Meter	15055
0.53mm	3.00	-20 to 260/280°C	30-Meter	15085

Rtx®-200 Structure



Visit www.restek.com for a complete listing of
Stabilwax® and **Rtx®-200** Columns.



New Stabilized PLOT Columns

Minimize Particle Release— Safe for Deans Switching Applications

NEW!**advanced
technology**

Restek has developed new procedures to manufacture PLOT columns with concentric stabilized adsorption layers. The new generation PLOT columns show a constant flow behavior and have significantly improved mechanical stability, resulting in easier operation, more reproducible retention times, and reduced particle release.

Rt®-Q-BOND Columns (fused silica PLOT) 100% divinylbenzene

ID	df (μm)	temp. limits	15-Meter	30-Meter
0.32mm	10	to 280/300°C	19743	19744
0.53mm	20	to 280/300°C	19741	19742

Rt®-QS-BOND Columns (fused silica PLOT)

porous divinyl benzene homopolymer

ID	df (μm)	temp. limits	15-Meter	30-Meter
0.32mm	10	to 250°C	19739	19740
0.53mm	20	to 250°C	19737	19738

Rt®-S-BOND Columns (fused silica PLOT) divinylbenzene 4-vinylpyridine

ID	df (μm)	temp. limits	15-Meter	30-Meter
0.32mm	10	to 250°C	19747	19748
0.53mm	20	to 250°C	19745	19746

Rt®-U-BOND Columns (fused silica PLOT)

divinylbenzene ethylene glycol/dimethylacrylate

ID	df (μm)	temp. limits	15-Meter	30-Meter
0.32mm	10	to 190°C	19751	19752
0.53mm	20	to 190°C	19749	19750

Rt®-Alumina BOND / Na₂SO₄ Columns (fused silica PLOT)

ID	df (μm)	temp. limits	30-Meter	50-Meter
0.32mm	5	to 200°C	19757	19758
0.53mm	10	to 200°C	19755	19756

Rt®-Msieve 5A Columns (fused silica PLOT)

ID	df (μm)	temp. limits	15-Meter	30-Meter
0.32mm	30	to 300°C	19720	19722
0.53mm	50	to 300°C	19721	19723

MXT®-Msieve 5A (Siltek®-treated stainless steel PLOT)

ID	df (μm)	temp. limits	30-Meter
0.53mm	50	to 300°C	79723



Got Sulfur?

**Sulfinert® treated
sample cylinders**

- Stable storage of sulfur and mercury at ppb levels
- Inert coating doesn't flake; more durable than Teflon®
- DOT rated to 1,800 or 5,000 psi.

www.restek.com/air



NEW!



Cylinder Holders, Wall Mounted

Description	qty.	cat.#
Single Cylinder Holder, Wall Mounted	ea.	21333
Double Cylinder Holder, Wall Mounted	ea.	23400
Triple Cylinder Holder, Wall Mounted	ea.	23401
Four Cylinder Holder, Wall Mounted	ea.	23402



NEW!

Swagelok® Stainless Steel & Brass Fitting kits

40% Savings by
purchasing kits!

Description	qty.	cat.#
Swagelok Brass Fitting Kit	kit	23141

Description	qty.	cat.#
Swagelok Stainless Steel Fitting Kit	kit	23197

Swagelok® Brass Fittings Kit

**Restek now carries a complete line of Swagelok® fittings!
For more information, visit us at www.restek.com**

Save a Life

(and your budget too!)

**Extend instrument lifetime and
reduce downtime for repairs.**

Routine maintenance essentials from Restek help you:

- Save money with LC and GC maintenance kits.
- Save time with innovative tools.
- Repair instruments no longer supported by the manufacturer.
- OEM performance at a fraction of the cost.



Handy kits and tools available at www.restek.com/savealife

contact your Restek representative today



Compare the Savings!

Inlet Liners for Agilent GCs

- Wide range of liners for Agilent GCs.
- Lower priced than other manufacturers.
- Performance equals or exceeds original manufacturer.

Description	Agilent cat.#		Restek cat.#	
Gooseneck Splitless (4mm) w/Wool (25-pk.)	5183-4694		22407	
4mm Split w/Wool (25-pk.)	5183-4692		20783	
4mm Split w/Wool (5-pk.)	5183-4691		20782	
Gooseneck Splitless (4mm) (25-pk.)	5183-4696		20800	
Gooseneck Splitless (4mm) (5-pk.)	5183-4695		20799	
Drilled Uniliner (hole near bottom) (ea.)	G1544-80730		20756	
Low Pressure Drop Precision Liner w/ Wool (5-pk.)	5183-4701		21033	
Gooseneck Splitless (4mm) w/Wool (5-pk.)	5183-4693		22406	
Double Gooseneck Splitless (4mm) (25-pk.)	5183-4706		20786	
Cup Splitter (5-pk.)	5183-4699		20710	

Check out
50% Savings
 on select liners for
 Agilent GCs at
www.restek.com/agilentliners

Custom Coatings for Ultimate Inertness!

Siltek® Deactivation for Liners

Our patented Siltek® deactivation process for liners produces a highly inert glass surface that features high temperature stability, extreme durability, and low bleed.

- Maximizes the inertness of the sample pathway with proprietary deactivation.
- Minimizes breakdown.
- Low bleed.
- Thermally stable.
- "Clean and green"—manufactured without the use of harmful organic solvents.
- Ideal for chlorinated pesticides.

For Siltek® deactivated inlet liners, add the corresponding suffix number to the liner catalog number.

qty.	Siltek Liner		Siltek Liner w/Wool		Siltek Liner w/CarboFrit	
each	-214.1	addl. cost	-213.1	addl. cost	-216.1	addl. cost
5-pk.	-214.5	addl. cost	-213.5	addl. cost	-216.5	addl. cost
25-pk.	-214.25	addl. cost	-213.25	addl. cost	-216.25	addl. cost

Call our technical service at 800-353-1688 ext. 4 for more information.

Restek can coat almost any chromatography part including:

- liners
- instrument parts
- sample loops
- fittings
- vials
- tubing



Lab Gas Issues?

ENSURE QUALITY AND MAINTAIN PRODUCTION WITH A RESTEK SOLUTION

GENERATION

Gas Generators

HYDROGEN



FID



ZERO AIR



NITROGEN



PURIFICATION

Super-Clean Gas Filter and Base Plate Kit

FREE Base Plate

with purchase of 2 triple or 2-fuel filters

Save
\$200

Description	qty.	cat.#
Restek Exclusive Triple Filter Kit: includes 2 triple filters & free base plate	kit	22859-409
Restek Exclusive Fuel Gas Filter Kit: includes 2 fuel gas filters & free base plate	kit	20529-409

Offer valid through Dec. 31, 2009.
Cannot be combined with other offers or discounts.





Good Deal!

Buy one 10-pk. of Inlet Seals Get one FREE!



Single-Column Installation (0.8mm Opening) Replacement Inlet Seal	Agilent part #	Restek cat.#	2nd 10-pk. Special Promo
Gold-Plated, 10-pk.	5188-5367	21318-409	FREE
Siltek Treated, 10-pk.		21320-409	FREE
Stainless Steel, 10-pk.	18740-20880	21316-409	FREE
0.25/0.32mm ID Dual-Column Installation (1.2mm Opening) Replacement Inlet Seal			
Gold-Plated, 10-pk.		21306-409	FREE
Siltek Treated, 10-pk.		21308-409	FREE
Stainless Steel, 10-pk.		20391-409	FREE
0.53mm ID Dual-Column Installation (1/16-inch Opening) Replacement Inlet Seal			
Stainless Steel, 10-pk.		20393-409	FREE
0.8mm ID Cross-Disk Inlet Seal			
Gold-Plated, 10-pk.	5182-9652	20467-409	FREE
Siltek Treated, 10-pk.		20474-409	FREE
1.2mm ID Cross-Disk Inlet Seal			
Gold-Plated, 10-pk.		21010-409	FREE
Siltek Treated, 10-pk.		21012-409	FREE

Offer ends **September 30, 2009**. Discount cannot be combined with any other offers.
Distributor participation may vary.

Great Deal!

Buy one 10-pk. of Dual Vespel® Ring Inlet Seals at 25% off list price & Receive a 2nd 10-pk. FREE!



**restek
exclusive!**

Dual Vespel® Ring Inlet Seals for Agilent GCs

0.8mm ID Dual Vespel Ring Inlet Seal	Restek cat.#	2nd 10-pk. Special Promo
Gold-Plated, 10-pk.	21241-409A	FREE
Siltek Treated, 10-pk.	21243-409A	FREE
1.2mm ID Dual Vespel Ring Inlet Seal		
Gold-Plated, 10-pk.	21247-409A	FREE
Siltek Treated, 10-pk.	21249-409A	FREE
0.8mm ID Dual Vespel Ring Cross-Disk Inlet Seal		
Gold-Plated, 10-pk.	22084-409A	FREE
Siltek Treated, 10-pk.	22086-409A	FREE

Dual Vespel Ring offer valid on Gold & Siltek inlet seals only.

Offer ends **September 30, 2009**. Discount cannot be combined with any other offers.

Distributor participation may vary.



Visit www.restek.com/dvr for all the features and benefits of this Restek innovation.

Your Column Installation Station

Make a clean cut and protect your column storage—it'll save you time and money!



20116



23015

Ceramic Scoring Wafer

Description	qty.	cat.#
Ceramic Scoring Wafers	5-pk.	20116

Scoring Wafer with Handle

Description	qty.	cat.#
Scoring Wafer with Handle	2-pk.	23015



22858



21044

Capillary Column Caps

- Attach to the column in seconds to form an airtight seal.
- Increase column lifetime—prevent moisture and air from entering the column during storage.

Description	qty.	cat.#
Capillary Column Caps (glass)	10-pk.	21044
Capillary Column Caps (silicone)	10-pk.	22858



New! Merlin Microseal Septa Kits NEW! for Varian and PerkinElmer GCs

Microseal Septa for PerkinElmer & Varian GCs	Merlin#	Similar to	cat.#
Varian 1078/1079 Kit (septum cap adapter, O-ring, nut, septum)	21-11w	—	22779
Varian 1177 Kit (septum cap adapter, O-ring, nut, spacer ring, septum)	22-11w	—	22780
PerkinElmer Kit (septum cap adapter, O-ring, nut, 2 septa)	51-12w	PE# N9303344	22781
Replacement Septum for Varian and PerkinElmer Kits (1 septum)	21-01w	PE# N9303345	22782

Did you know that Restek carries a complete line of ETP Electron Multipliers for Mass Spectrometry?



Description	qty.	cat.#
Electron Multipliers for Agilent GC/MS and LC/MS		
For Agilent 5970 GC/MS	ea.	23072
For Agilent 5971, 5972, GC GC/MS	ea.	23073
For Agilent 5973 & 5975 GC/MS (includes mount for initial installation)	ea.	23074
For Agilent 5973 & 5975 GC/MS and LC/MSD (Replacement Multiplier)	ea.	23075
For Agilent LC/MSD (includes mount for initial installation)	ea.	23076
Electron Multiplier for Applied Biosystems (Sciex)		
For API 300, 3000 & 4000 Applied Biosystems	ea.	23077
Electron Multiplier for Thermo Finnigan GC/MS		
For Thermo TRACE DQ, DSQII, and Polaris-Q GC/MS	ea.	23081

Other electron multipliers are available upon request. If you don't see your instrument listed, give us a call at 814-353-1300 ext. 3



QuEChERS

Quick, Easy, Cheap, Effective, Rugged and Safe!
dSPE Tubes for pesticide residue analysis from food products

- Fast, simple sample extraction and cleanup using dSPE.
- Fourfold increases in sample throughput.
- Fourfold decreases in material cost.
- Convenient, ready to use centrifuge tubes with ultra pure, preweighed adsorbent mixes.



FREE

Sample Packs
Available!

**For Greater
Convenience**

NEW Standards available!

Complete Product Offering

**to support AOAC (2007.01)
& European (EN15662)
QuEChERS Methods**

- Comprehensive portfolio of QuEChERS tubes
- Method-specific internal and QC standards
- Rxi®-5Sil MS Capillary column **on sale now—25% off!**
- Ultra Aqueous C18 LC Column

View our QuEChERS On-Demand Webinar!



www.restek.com/quenchers

Extraction Cells for ASE® 200



Buy 5 cells or cap pairs & get the 6th one free!

- Meets original equipment manufacturer's performance.
- Choose original equipment-equivalent stainless steel, or Siltek® deactivation for improved inertness and extended lifetime. (Discount applies to stainless steel only.)

Save now until June 30, 2009 on these popular sizes:

Description	qty.	Stainless Steel cat.#
Extraction Cell Body for ASE 200, 22mL	six for the price of 5	26098-600
Extraction Cell Body for ASE 200, 33mL	six for the price of 5	26116-600
Replacement Extraction Cell End Caps for ASE 200	six 2-pks. for the price of 5	26096-600

Syringe Filters with Luer Lock Inlet



100 pack

Cut costs, not corners!

Top quality filters at lower costs than even the discounted prices of other manufacturers!

Compare our quality with your current brand using our **FREE no risk sample packs. Call today!**

PTFE Purple Luer Lock Inlet 13mm, 0.22µm, 100pk

Request lit. cat# GNFF1072 for a membrane selection guide and 12 reasons to buy Restek syringe filters.

Visit www.restek.com/filters for complete product listing.



Free Data Packs available online!

www.restek.com/datapacks



GET YOUR MIX

Time-Saving MegaMix® 
Environmental Reference Mixes.

- Largest number of target analytes in one mix, formulated for maximum stability.
- Available for US EPA methods 8260, 8270, 502.2, 524.2, 525.2, 624, 625, SOM01.1, OLC 03.2, OLM 04.2, Skinner List volatiles, Skinner List semivolatiles.

Save time, save effort, minimize potential for preparation problems—use MegaMix® reference mixes, only from Restek or authorized distributors.



www.restek.com/megamix



Big Savings with Restek Mini-Cans and sampling kits for vapor intrusion



- Smaller cans are less expensive to buy and cheaper to ship for analysis.
- Quick-Connect Fittings provide fast, easy can-to-kit connections in the field—no tools required!

Miniature Air Sampling Canisters with Quick-Connect Fitting

Description	qty.	400cc		1,000cc
		cat.#		cat.#
Siltek Treated Miniature Canister with Siltek Treated Quick-Connect Stem Fitting	ea.	24190		24196

Quick-Connect Fittings for Miniature Air Sampling Canisters

Description	qty.	cat.#
Quick-Connect Body Fitting	ea.	24187

Note: Quick-connect body fitting (cat.# 24187) must be ordered separately to sample with quick-connect stem fitting.

Miniature Air Sampling Kits

NEW!

Canister Volume*/Sampling Time	Flow (sccm)	Orifice size	Siltek Treated Sampling Kits
400cc	1 Liter		
8 hour	24 hour	0.5–2	0.0008" 26253
1 hour	2 hour	4–8	0.0016" 26257

*Additional orifice sizes available.



get mini!

For more information on Restek's Mini-Cans and sampling kits, visit www.restek.com/air

Fill Tedlar® Bags without Pump Contamination with Vacuum Bag Sampler

Buy sampler & Get pack of 10 bags FREE!

Save Now!

- Cost-effective alternative for whole air sampling.
- Light weight and easy to use.
- Vacuum sampler fills bag quickly, without sample passing through pump.
- Large selection of Tedlar® bag sizes available, as well as stainless steel valve option.

Description	qty.	cat.#
Bag Sampler/Tedlar Bag		
Includes Vacuum Bag Sampler Model 1062 and 10pk 1L Tedlar Bags w/PP valve (22118 + 22050)	kit	22118-406.1
Bag Sampler/Tedlar Bag		
Includes Vacuum Bag Sampler Model 1062 and 10pk 10L Tedlar Bags w/PP valve (22118 + 22053)	kit	22118-406.2



Offer valid through June 30, 2009.



NEW!

Ultra II LC Columns

Restek Manufactured, High Purity, Type B Silica for Column-to-Column Reproducibility for ANY LC System

Available phases:

- C18
- Silica
- Biphenyl
- PFP Propyl
- Aromax
- Aqueous C18

Available Particle Sizes:

- 1.9µm for UHPLC
- 2.2µm for UFLC and RRLC
- 3µm and 5µm for HPLC

More phases coming soon!

Save on Preventive Maintenance Kits & Parts for Agilent & Waters Systems

Description for Agilent LCs	Model #	Similar to Agilent		
		part #	qty.	cat.#
Preventive Maintenance Kit	1050	01078-68721	kit	25259
Preventive Maintenance Kitt cap	1050	—	kit	25259A
Autosampler Preventive Maintenance Kit	1100	G1313-68709	kit	25271
Pump Maintenance Kit	1050, 1100	G1311-68710	kit	25270

Description for Waters LCs	Model #	Similar to Waters		
		part #	qty.	cat.#
Preventive Maintenance Kit	Alliance 2690, 2695	WAT270944	kit	25143
Pump Preventive Maintenance Kit	600 Pump	WAT052675	kit	25144
Autosampler Preventive Maintenance Kit	717 Autosampler	WAT052669	kit	25145
Pump Preventive Maintenance Kit	616 Pump	WAT052672	kit	25146
Pump Preventive Maintenance Kit	515 Pump M6KA, 501, 510, 515, 590,	WAT052587	kit	26519
PerformancePlus Check Valve Cartridge	600E	700000254	2-pk.	25370




Replacement parts are also available for PerkinElmer, Shimadzu, Thermo Scientific, Beckman & Hitachi HPLC Systems at very affordable prices.



Don't burn up your budget!

Save now on selected Detector Lamps.



Description	Similar to part #	qty.	cat.#
For Agilent HPLC Systems			
Detector Lamp, 1090 DA, 1050 VW/DA/MWD, for Agilent Model #s 1090, 1050	79883-60002	1-pk. 2-pk. 3-pk.	25260-712.10 25260-712.15 25260-712.20
Lamp, VWD G1314A, for Agilent Model #s 1100, 1200	G1314-60100	1-pk. 2-pk. 3-pk.	25262-712.10 25262-712.15 25262-712.20
8453 Deuterium Lamp, for Agilent Model # 8453	2140-0605	1-pk. 2-pk. 3-pk.	25263-712.10 25263-712.15 25263-712.20
Lamp, DAD Long Life Deuterium (2000 hours), for Agilent Model #s 1100	5181-1530	1-pk. 2-pk. 3-pk.	25399-712.10 25399-712.15 25399-712.20
For Shimadzu HPLC Systems			
Deuterium Lamp, for Shimadzu Model #s SPD-10A, 10AV, M20A	228-34016-02	1-pk. 2-pk. 3-pk.	25284-712.10 25284-712.15 25284-712.20

Universal 10-32 PEEK™ Column Connectors and Plugs

Universal PEEK™ Connectors allow easy installation of all 1/16-inch tubing, including stainless steel.



Description	qty.	cat.#
PEEK Column Connector (beige, round body)	10-pk.	25015
PEEK Column Plug (black)	10-pk.	25016
PEEK Fingertight Fittings (blue, flat-sided)	10-pk.	25324

Low-Pressure Slip-On Inlet Filter for Mobile Phase Reservoir



This slip-on 10µm porosity filter easily attaches to the pump inlet line, without the use of tools. The universal SST tip with Tefzel® collar accommodates standard Teflon® tubing inner diameters. 1/8" OD fits Altex, ISCO, LDC, Varian, Waters, PerkinElmer, and other pumps.

Description	qty.	cat.#
Slip-on Inlet Filter	ea.	25008

Why allow a tiny leak to become a GIANT problem?

Restek Electronic Leak Detector

Protect your data and
analytical column!

NEW!

- A sleek, new ergonomic, hand-held design.
- Longer lasting battery, up to 6 hours of continuous use.
- Automatic shut-off capabilities.
- A convenient carrying and storage case.
- A universal charger set (US, European, UK, and Australian plugs included).

Leak Detector with Universal Charger Set
(US, UK, European, Australian)

Cat.# 22839

For more information, visit
www.restek.com/leakdetector

 RESTEK

Lit. Cat.# GNES1147-INT
2009 v.4

© 2009 Restek Corporation.

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823 • 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: +33 (0)1 60 78 32 10 • fax: +33 (0)1 60 78 70 90 • e-mail: restek@restekfrance.fr

Restek GmbH • phone: +49 (0)6172 2797 0 • fax: +49 (0)6172 2797 77 • e-mail: info@restekgmbh.de

Restek Ireland • phone: +44 (0)2890 814576 • fax: +44 (0)2890 814576 • e-mail: restek@restek.com

Restek Japan • phone: +81 (3)6459 0025 • fax: +81 (3)6459 0025 • e-mail: ryosei.kanaguchi@restek.com

Thames Restek U.K. LTD • phone: +44 (0)1494 563377 • fax: +44 (0)1494 564990 • e-mail: sales@thamesrestek.co.uk

ISO 9001:2000
cert.# FM80397